

TruBio Product Family

Features	TruBio PLSS or RDPD	TruBio SCADA
Basic modules and graphics	Yes	SCADA loop blocks
TruConfig	Yes	Yes
<ul style="list-style-type: none"> ■ ConfigBuilder ■ Vessel Save&Load ■ Signal Conditioning ■ Offline Data Entry 	Data entry with auto standardize	Data entry only
TruCal	Yes	Yes
<ul style="list-style-type: none"> ■ Text Prompts Guide ■ Electronic Report 		
TruHealth	Yes	If data is available
<ul style="list-style-type: none"> ■ Secondary Measurement Signals 		
Simulate	Yes	N/A
TruConnect	Yes	Yes
<ul style="list-style-type: none"> ■ Functions can be restricted according to the TruBio licenses 		
Plug-In		
<ul style="list-style-type: none"> ■ DataGrabber ■ DataConnector ■ OSI Pi 	Nova BioProfile 400 Nova FLEX	



Finesse Solutions, LLC
3350 Scott Blvd, Bldg 1
Santa Clara, CA 95054

Finesse Solutions AG
Via Sogn Gieri 27a
CH-7402 Bonaduz
Switzerland

877-204-8644
www.finesse.com